




SYNERGY FLEX

REVEAL BEAUTIFUL SKIN

SYNERGY FLEX
cacibeauty.com

CACI



WORLD LEADERS
IN AESTHETIC
TECHNOLOGY

ABOUT CACI INTERNATIONAL

CACI International has revolutionised the beauty industry since its beginning in 1992, providing cutting edge technology and maintaining an incomparable reputation for excellence and innovation. We've paved the way in anti-ageing, skin rejuvenation and skin solutions for over three decades.

As the UK's leading equipment manufacturer and supplier of non-invasive aesthetic treatment systems, we have supplied to over 10,000 of the world's leading hotels, spas and aesthetic clinics in over 100 countries.

We use the latest technological and scientific advancements to achieve the very best treatment results that can be tailored for men and women of all ages and skin types.

Our treatments are often described as the 'red-carpet beauty secret' and have attracted an impressive celebrity following including royalty, film stars, supermodels, editors and make-up artists from around the globe.

THE CACI DIFFERENCE

Whilst CACI systems incorporate multiple technologies, at the heart of all our systems is our core microcurrent technology. What sets our microcurrent technology apart from other microcurrent machines is our use of ultra-low frequencies in combination with skin resistance technology and a patented digital waveform.

This unique digital waveform optimises results and comfort by gradually increasing in intensity to overcome skin resistance. It then rises to a peak to deliver the full power of the microcurrent energy deep into the muscle. This means that CACI systems not only achieve better results but also longer lasting results.

With the introduction of the **Synergy Flex system**, we have taken Microcurrent technology to the next level.

The multi-functional Synergy Flex is the latest innovation from CACI, it is our most advanced facial system, clinically proven to **reduce wrinkle depth by up to 75% and improve skin elasticity by up to 88%.**

S.P.E.D® DUAL ACTION TECHNOLOGY


The system offers a range of exclusive technologies, including our unique **S.P.E.D® Dual Action Technology** that provides simultaneous skin rejuvenation and facial toning. Synergy Flex is clinically proven* to provide even faster, better, and longer-lasting results than our earlier CACI facial systems.

CACI is a brand synonymous with delivering effective treatment results without the need for surgery. When you invest in a CACI, you join a winning team and share in the success of the CACI brand. A CACI system will attract new clients, drive treatment revenues, and help your business grow.

We also support your business by providing a wide range of promotional and marketing materials, including a marketing guide with tips on how to promote your business, posters, treatment brochures, window stickers, and a variety of downloadable online resources.

Our dedicated team of business consultants is ready to advise and help your business succeed.

At CACI, we are committed to being your partner.



THE CACI
DIFFERENCE

CUSTOMISABLE SYSTEM. FLEXIBLE TECHNOLOGIES.

The Synergy Flex system offers a range of exclusive technologies for skin rejuvenation and facial toning. The Synergy Flex comes equipped with all CACI advanced attachments, elevating your business and therapists to Advanced Clinic/Therapist status. As an advanced clinic you will receive enhanced website listings, advanced clinic recognition, and access to premium point-of-display materials.

VERSATILE AND FUTURE-PROOF DESIGN

Discover the new look and enhanced versatility of the Synergy Flex, future-proofed with its independent Flex Tec and Microcurrent Innovation ports. The Flex Tec port provides

a dedicated 12V power source for future technology handsets, ensuring your Synergy Flex system stays current as new innovations emerge.

The Microcurrent Innovation port features a specialised microcurrent power source, driving future microcurrent advancements. The customisable screen allows users to select any frequency, waveform, or power level, and save these settings as custom presets, offering unparalleled flexibility in microcurrent treatments.

Additionally, the system includes three independent microcurrent channels, allowing for simultaneous treatments. This means, for example, a client's face and hands can be treated at once, boosting business revenue without extending treatment time.

With Synergy Flex, your system remains cutting-edge, seamlessly integrating new technologies as they become available.

WHY CHOOSE SYNERGY FLEX?

Innovative Microcurrent Technology: Achieve superior and longer-lasting results with our unique digital waveform and ultra-low frequency combination.

S.P.E.D® Dual Action Technology: Deliver the power of LED light stimulation at the same time as using Microcurrent. This dual action technology has been clinically proven to deliver faster, more effective, and longer lasting results than standard microcurrent.

Exclusive Skin Rejuvenation Technologies: Enhance your offerings with advanced technologies specifically designed for optimal skin rejuvenation.

Future-Proof Flexibility: Stay ahead with independent Flex Tec and Microcurrent Innovation ports for seamless integration of new technologies.

Comprehensive Treatment Programs: Offer a wide variety of treatments tailored to different skin concerns and types.







MICROCURRENT DUAL PROBES WITH S.P.E.D® TECHNOLOGY

The dual-tipped Microcurrent Probes utilise CACI's pioneering S.P.E.D® Dual Action Technology, which allows skincare specialists to simultaneously utilise LED light therapy and microcurrent. This innovative dual-action approach enhances the effectiveness of treatments, providing clients with superior results.

The Microcurrent Dual Probes are designed to deliver targeted microcurrent therapy while emitting LED light, maximising the benefits of both technologies. This combination stimulates tissue regeneration, boosts collagen production, and tightens sagging muscles, all while smoothing the appearance of fine lines and wrinkles. The synergistic effects of LED light and microcurrent energy work together to provide clients with a comprehensive solution for achieving firmer, smoother, and more youthful-looking skin.



MICROCURRENT QUAD PROBES WITH S.P.E.D® TECHNOLOGY

The Microcurrent Quad Probes are engineered to grip a larger surface area and effectively target the jaw muscles, particularly within the Jowl Lift treatment. This specialised treatment addresses muscle laxity around the jawline. By focusing on this area, the treatment lifts and firms the muscles, redefining facial contours to give a firmer, more toned appearance.

The Quad Probes, equipped with S.P.E.D® Dual Action Technology, enhance the treatment by combining microcurrent therapy with LED light stimulation. This dual-action approach accelerates the lifting and toning process, providing faster and more effective results. The LED light not only boosts the skin's natural healing processes but also promotes collagen production, resulting in improved skin elasticity and a more youthful look.



MICROCURRENT HYDRATONE ROLLERS WITH S.P.E.D® TECHNOLOGY

Hydratone S.P.E.D® Dual Action Technology Rollers are designed to deliver unparalleled skin hydration and rejuvenation. The rollers gently massage the face over the electrically conductive CACI Hydro Mask or Purifying Silver Mask, transmitting the rejuvenating energy of microcurrent deep into the skin to restore muscle tone and firmness. The red LED light is primarily used for its anti-aging properties. It stimulates collagen production, which is essential for maintaining skin elasticity and firmness. This rejuvenating effect helps to reduce the appearance of fine lines and wrinkles, resulting in a smoother, more youthful complexion. The blue LED light has a calming effect on the skin, making it ideal for treating inflamed or irritated areas. It helps to soothe the skin, reduce redness, and promote a clearer complexion.



MICROCURRENT EYE ROLLERS

The Eye Revive Rollers are designed for treating the delicate eye area. Eye Revive Gel is applied directly into the roller hand piece, these rollers enhance product absorption into the skin. Utilising microcurrent technology, the Eye Revive Rollers gently tighten and tone sagging muscles around the eyes. Paired with the deeply nourishing CACI Eye Revive Mask, this treatment lifts hooded eyelids, diminishes puffiness and dark circles, and visibly softens the appearance of fine lines and wrinkles around the eyes.

The Eye Revive Rollers provide a treatment solution for rejuvenating the delicate eye area, offering clients a rejuvenated and refreshed look with each treatment session.



MICROCURRENT ELECTRO GLOVES

The Electro Gloves are made from a unique silver fabric that is electrically conductive.

During the Micro-Touch Facial, the Electro Gloves are used to deliver microcurrent electrical impulses that help to improve skin firmness and boost circulation. This innovative approach allows for precise targeting of facial muscles and skin tissues, enhancing the overall firmness and elasticity of the skin.

Within the Hand Rejuvenation Treatment, the Electro Gloves are fitted over a Rejuvenating Hand Mask which intensively moisturises the skin. The combination of the Electro Gloves and the hand mask ensures that microcurrent energy is effectively delivered to the deeper layers of the skin. This treatment leaves the clients' hands hydrated, with a softer, smoother, and more youthful appearance.



MICROCURRENT ENERGISING MASK

The Microcurrent Energising Mask is an innovative sheet mask featuring breakthrough silver ion conductive print technology. The mask is printed with an electronic circuit that when activated by the Synergy Flex delivers microcurrent energy to rejuvenate the skin. The mask is attached to the Synergy Flex using a cable with magnetic connectors.

Each mask is single use and can be used as a stand-alone 15-minute skin treatment or seamlessly integrated into a CACI facial treatment.



MICROCURRENT GUA SHA

The Microcurrent Gua Sha is an optional microcurrent innovation attachment exclusive to Synergy FLEX. It merges the benefits of microcurrent therapy with traditional Gua Sha techniques.

This advanced skincare tool sculpts and contours the face and body, reduces fine lines, and promotes lymphatic drainage for a healthier, more radiant glow. It effectively relieves tension in muscles, enhances skin texture and tone, and offers a deeply relaxing experience while delivering the powerful benefits of microcurrent technology.

ENHANCING MICROCURRENT TREATMENTS



CONDUCTIVITY BOOSTERS WITH SILVER FIBRE

The Conductivity Boosters, made from highly conductive silver fibres, are inserted into the Microcurrent Probes to enhance the flow of microcurrent, maximising treatment results.

The Boosters offer superior skin grip and a powerful lift, while their natural anti-bacterial properties help clarify the skin and reduce breakouts. The serum inside the boosters is infused with electrolytes, hyaluronic acid, niacinamide, and peptides which boosts conductivity, hydrates, and rejuvenates the skin. For optimal results, pair with our Super Charge Energisers during microcurrent treatments.



SUPER CHARGE ENERGISERS

The Super Charge Energisers are sachets of dissolvable powder, specifically formulated to enhance microcurrent treatments and improve microcurrent conductivity. They have been designed to be used with CACI Microcurrent Probes, Microcurrent Energising Masks and Electro Gloves.

The Super Charge Energisers boost the electrical conductivity of water by a staggering 196% (results based on one sachet mixed with 120ml of water). The Energisers also include Vitamin C, Marine Collagen and Niacinamide to leave your clients' skin feeling super soft and silky smooth.



ORBITAL ABRADER

The Orbital Abrader uses an exfoliation technology with a hand piece featuring an abrasive tip that rotates at varying speeds on an elliptical axis. This orbital rotation performs controlled skin exfoliation, removing imperfections and impurities. Single-use, disposable abrasive tips coated with silicon carbide allow for either superficial or deep exfoliation of the epidermis, resulting in a clearer, smoother, and more youthful complexion. Precision control over the depth of abrasion lets you tailor each treatment to the client's needs.

This crystal- and vacuum-free method of microdermabrasion avoids the common blockage, maintenance, and downtime issues associated with crystal-based systems. Additionally, a massage tip can be attached to the hand piece for lymphatic drainage treatments. The Orbital Abrader also delivers simultaneous LED light during exfoliation to promote healing and skin rejuvenation.



ULTRASONIC ACTUATOR

The Ultrasonic Actuator hand piece contains ceramic/quartz transducers that when excited by an electrical current emit rapid vibrations, 27,000 per second. The vibrations propel impurities to the skin surface and cause dead skin cells to become dislodged. The hand piece is then used to deeply cleanse and gently exfoliate the epidermis. The Ultrasonic Actuator provides a far gentler and less aggressive method of skin exfoliation than traditional microdermabrasion.

The hand piece can also be set to a specific frequency to enhance the absorption of skin care products using sonophoresis (ultrasound waves), to push the products deeper into the skin.



WRINKLE COMB

The Wrinkle Comb is used to effectively and painlessly plump out and soften deep lines and wrinkles, thereby providing a needle-free alternative to dermal fillers.

The unique handset transmits high frequency electrical stimulation simultaneously with LED light, which stimulates an increase in the production of elastin and collagen. Red LED light encourages skin rejuvenation and promotes tissue repair. Blue LED Light has an antibacterial and calming effect on inflamed or irritated skin. An alternating pulsed LED mode is particularly effective at treating acne-prone or problematic skin.



INDEPENDENT CLINICAL STUDY

CACI Synergy Flex, introduces pioneering S.P.E.D® Dual Action Technology that delivers simultaneous LED skin rejuvenation and microcurrent facial toning. An independent clinical study conducted by Cutest Systems Ltd, compared CACI Synergy Flex S.P.E.D® enhanced microcurrent with standard microcurrent facial toning. 25 female participants aged between 35-60 took part in the study, 10 treatments were administered with measurements and assessments taken after treatment 1, 5, 10 and 2 weeks following the final treatment.

RESULTS

STUDY PARTICIPANT QUESTIONNAIRE

84%

Agreed skin felt firmer

92%

Agreed skin felt smoother

80%

Agreed fine lines reduced

88%

Agreed saw a visible improvement

SKIN ELASTICITY IMPROVEMENT

Up to

66%

Standard Microcurrent

VS

Up to

88%

CACI S.P.E.D®

After 10 treatments

WRINKLE DEPTH REDUCTION

47%

Standard Microcurrent

VS

70%

CACI S.P.E.D®

After 5 treatments

SUMMARY


CACI S.P.E.D® provides a greater reduction in red and brown spots than standard microcurrent.

CACI S.P.E.D® reduces the appearance of blemishes on the skin more effectively than standard microcurrent.


CACI S.P.E.D® results in a faster improvement in skin texture than standard microcurrent.

CACI S.P.E.D® produces faster, better and longer lasting results in the improvement of lines and wrinkles compared to standard microcurrent.


CACI S.P.E.D® produces faster, better and longer lasting results in improving skin elasticity compared to standard microcurrent.




BEFORE



BEFORE



BEFORE



AFTER



AFTER



AFTER

INCREASE YOUR REVENUE

WITH TREATMENTS THAT APPEAL TO ALL CLIENTS

Choosing a CACI Synergy Flex system gives you the comfort and peace of mind knowing that you are purchasing direct from the UK's largest aesthetic equipment manufacturer and from a brand that you can trust to deliver quality, reliability and performance.

CACI Synergy Flex offers a wide range of treatment options that will appeal to men and women of all ages and skin types.



Offering CACI treatments within your business will enable you to attract new clients, increase client retention and significantly boost your treatment revenues.

With a variety of attractive payment schemes available you don't even need to make a capital outlay and can start generating profit as soon as your CACI Synergy Flex system arrives.

Our payment schemes are so readily affordable that you only need to perform just one CACI client treatment course per month to be in profit.

A CACI Synergy Flex system really is a must have for any beauty business serious about generating substantial additional income.

SHARE IN THE SUCCESS OF THE CACI BRAND NAME AND INCREASE YOUR REVENUE



Based on a 7 day week and £135 per CACI Synergy Flex 90 minute Treatment

“ I INTRODUCED CACI INTO MY SALON 23 YEARS AGO AND HAVE NEVER LOOKED BACK. WE HAVE A VERY LARGE CACI CLIENTELE, AND HAVE UPDATED OUR CACI SYSTEMS OVER THE YEARS. OUR CLIENTS ARE DELIGHTED WITH THE RESULTS WE CAN ACHIEVE.”

Ashleigh Bland, Ponteland Hair & Beauty, Newcastle upon Tyne.



SYNERGY FLEX

SIGNATURE TREATMENTS

CACI Signature Non-Surgical Facial Toning **(60 minutes) RTP £90.00**

This non-invasive facial features S.P.E.D® Dual Action Technology that provides simultaneous skin rejuvenation and facial toning. This treatment uses microcurrent impulses, to lift and tone and LED light therapy to improve skin elasticity and reduce the appearance of fine lines and wrinkles. Firms the face and gives the skin a more youthful, fresh appearance and glow.

CACI Synergy Advanced Non-Surgical Facial Toning **(90 minutes) RTP £135.00**

Combines CACI's Signature Non-Surgical Facial Toning with advanced skin exfoliation techniques to revitalise the skin, resulting in a brighter and smoother complexion. A combination of advanced technologies are used to plump fine lines and wrinkles, instantly lift and tone the face and firm the neck.

CACI Synergy Purifying **(45 Minutes) RTP £60.00**

This treatment tackles problematic skin and leaves the complexion looking smooth and radiantly healthy. A combination of red and blue LED light therapy, deep cleansing and skin exfoliation techniques are used to effectively brighten and even the skin tone.

CACI Gentleman's Facial **(45 minutes) RTP £65.00**

Using a combination of advanced technologies, this treatment is designed to specifically address male skincare concerns. Deep cleansing and skin exfoliation techniques tackle razor bumps and in-grown hairs, electrical impulses tighten, firm and redefine the jawline; the CACI Hydro Mask and LED light therapy are then used to hydrate, calm and smooth the skin, leaving the complexion refreshed, revitalised and more youthful looking.

CACI Skin Calm **(30 minutes) RTP £40.00**

Targets skin prone to breakouts and irritation. This treatment uses LED light therapy that is renowned for its anti-bacterial and anti-inflammatory healing properties, combined with a gel mask to soothe, calm, reduce redness and rehydrate the skin.

ADVANCED TREATMENTS

CACI Flex Facial* **(60 minutes) RTP £120.00**

The Flex Facial is a transformative treatment combining CACI's pioneering microcurrent technologies. Targeting muscle laxity around the jawline with S.P.E.D® Dual Action Technology, this cutting-edge treatment harnesses the power of LED light therapy alongside microcurrent. It's followed by the innovative CACI Microcurrent Gua Sha, designed to promote lymphatic drainage, boost circulation, and reduce puffiness. The result is lifted, firm, and beautifully contoured facial features. The experience concludes with the Microcurrent Energising Mask, renewing and repairing the skin while stimulating collagen and elastin for a radiant, glowing complexion.

CACI Microcurrent Gua Sha* **(45 minutes) RTP £70.00**

Combining CACI Microcurrent with traditional Gua Sha technique, where high-tech innovation meets holistic therapy. This treatment rejuvenates your skin by promoting lymphatic drainage, boosting circulation, and reducing puffiness, whilst using the CACI's microcurrent to contour and sculpt the facial features.

CACI Jowl Lift **(15 minutes) RTP £35.00**

Target muscle laxity around the jawline using S.P.E.D® Dual Action Technology, which harnesses the power of LED light therapy at the same time as using CACI microcurrent. This treatment will lift and firm the muscles and redefine facial contours to give a firmer, more toned appearance. Ideal for women wanting to improve the appearance of sagging jowls and for men who desire a more chiselled jawline.

CACI Eye Revive **(30 minutes) RTP £45.00**

Lift and firm muscles around the eye area, particularly where you may tend to lose definition as you age and develop 'hooded eye lids'. This treatment uses serum filled microcurrent rollers together with the soothing and calming Eye Revive Mask to reduce puffiness and dark circles, and soften the appearance of fine lines and wrinkles around the eyes.

CACI Hand Rejuvenation **(30 minutes) RTP £40.00**

The anti-ageing and nourishing hand treatment uses conductive microcurrent gloves with CACI Rejuvenating Hand Mask to reduce fine lines and intensively moisturise, leaving the hands hydrated with a softer more youthful appearance.

CACI Micro-Touch Facial **(45 minutes) RTP £45.00**

The Micro-Touch Facial combines the muscle toning and lifting benefits of CACI's iconic microcurrent technology with hands-on treatment therapy using the Electro Gloves. The treatment focuses on muscle lift, tension release, deep tissue relaxation, and lymphatic drainage for more sculpted, toned, enhanced facial definition.

EXPRESS TREATMENTS

CACI Hydratone **(15 minutes) RTP £45.00**

This treatment is designed to deeply nourish, hydrate and firm the skin. Facial rollers emitting simultaneous microcurrent impulses and LED light therapy are gently massaged over a hydrating gel facial mask. The rejuvenating energy transmitted through the rollers provides immediate and intensive hydration and helps to restore skin firmness, soften fine lines and rejuvenate the skin. Ideal for the treatment of dehydrated or sun damaged skin.

CACI Wrinkle Revolution **(30 minutes) RTP £40.00**

Targets deep lines and wrinkles, such as crow's feet and laughter lines. This treatment is a non-invasive, needle free alternative to collagen injections and dermal fillers that will instantly plump fine lines and wrinkles from the very first treatment.

CACI Lip Booster **(15 minutes) RTP £25.00**

Enhance your pout. Designed to give fuller, plumper lips and soften top lip lines. This quick fix treatment will have your clients' pout ready in no time.

CACI Skin Confidence **(30 minutes) RTP £40.00**

Designed to soften and refine the appearance of stretch marks and scarring, this treatment uses advanced skin exfoliation techniques together with the regenerating effects of microcurrent stimulation and LED light therapy to trigger skin tissue repair and healing.

CACI Re-Energise & Renew **(15 minutes) RTP £30.00**

This facial will revitalise and re-energise your skin. The treatment infuses your skin with microcurrent energy using the Microcurrent Energising Mask which boosts blood flow and stimulates tissue repair. It helps reduce blemishes, pigmentation, and smooths out fine lines and wrinkles, while promoting collagen and elastin production for firmer, more radiant skin. Perfect for rejuvenating your complexion, this facial leaves your skin looking renewed and glowing.

Please note: Express treatments can be performed as either stand alone or as add ons to any of our other treatments.

Please note: Advanced & Express Treatments can be performed as either stand alone or as add ons to any of our other treatments.

*These treatments can only be performed using the optional advanced attachments.



TARGETED RETAIL SKINCARE



TECHNICAL SPECIFICATIONS

SYSTEM

Dimensions: **W38 x D40 x H24cm**
 Weight: **4.5kg**
 Touch Screen: **11.6in**
 Supply Voltage: **90-250 Volts**
 Output Voltage: **30 Volts max (loaded)**
 Consumption: **110V pk-pk max**
 Safety class: **1 type BF in accordance with IEC 601-2-10**

MICROCURRENT

Channels: **3 independent**
 Polarity: **Positive, negative & alternating**
 Range: **100uA to 640 uA**
 Frequency: **0.3Hz-300Hz**
 Waveshapes: **12**
 LED BLUE: **470nm - 60mW**
 LED RED: **620nm - 60mW**

ULTRASONIC ACTUATOR

Frequency: **22,000Hz - 28,000Hz**

FLEX TEC

Output: **12V DC**

WRINKLE COMB

Output: **0.2mA - 1.6mA**
 Output Frequency: **6000Hz**
 LEDs: **5x Blue 470nm-40mW**
 LEDs: **5x Red 620nm-30mW**

ORBITAL ABRASION

Rotation Speed: **0-4000 rpm**
 Motor Voltage: **6-15 Volts**
 Amperage: **0.1-0.5 Amps**
 LEDs: **6x Blue 470nm-40mW**
 LEDs: **6x Red 620nm-30mW**

Pictured Synergy Flex with optional trolley.



CACI International, 514 Centennial Avenue,
 Centennial Park, Elstree, Hertfordshire, WD6 3FG, UK

T +44 (0)20 8731 5678
 E info@caci-international.co.uk



CACI